

Technical Datasheet

		EN ISO 13485 certified
Article - No:	10-131-0000	Revision: 2
Remarks		single use only

Specification

Product	Lid PS for micro test plate 96 wells, highly transparent, np pcr ready
Dimensions	length: 126,90 mm (+/- 0,2 mm)
	width: 85,00 mm (+/- 0,2 mm)
	height: 9,70 mm (+/- 0,2 mm)
Colour	highly transparent
Material Characteristics	polystyrene (PS)
Temperature range	-20°C - 60°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 5
	pcs per carton: 100
	pcs per pallet: 6.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

