

Technical Datasheet

EN ISO 13485 certified	
Article - No:	11-022-0000
Remarks	Revision: 2 single use only

Specification

Product	Mailing container PP, 40ml, Ø30x85 mm, with absorbing inlay and screw thread, without cap, natural, autocl. up to 121°C
Dimensions	diameter: 30,00 mm (+/- 0,2 mm)
	height: 85,00 mm (+/- 0,2 mm)
Colour	natural
Nominal Capacity	40ml
Material Characteristics	container: polypropylene (PP)
	inlay: cellulose
Temperature range	-196°C - 121°C
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 500
	pcs per carton: 500
	pcs per pallet: 12.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

