

Technical Datasheet

EN ISO 13485 certified	
Article - No:	11-162-0000
Revision:	2
Remarks	single use only

Specification

Product	Mailing container PP, 60ml, Ø30x126 mm, with absorbing inlay and screw thread, without cap, natural, autocl. up to 121°C
Dimensions	diameter: 30,00 mm (+/- 0,2 mm)
	height: 126,00 mm (+/- 0,2 mm)
Colour	natural
Nominal Capacity	60ml
Material Characteristics	polypropylene (PP)
Temperature range	-196°C - 121°C
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 300
	pcs per carton: 300
	pcs per pallet: 7.200
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

nerbe plus GmbH & Co. KG
 Porschestr. 25
 D-21423 Winsen/Luhe
 Tel: +49 (0)4171 7849 0
 Fax: +49 (0)4171 6 45 45

info@nerbe-plus.de
 www.nerbe-plus.de
 HRA Lüneburg 203240
 USt.-ID: DE 32583 0652

Pers. haftende Gesellschafterin:
 NERBE R&D GmbH, Winsen/Luhe
 GF: Philipp J.F. Nerbe
 HRB Lüneburg 208826

Deutsche Bank AG
 Kto. 0204 41200
 IBAN: DE41 2407 0075 0020
 4412 00
 SWIFT / BIC DEUTDE2240

Hamburger Sparkasse
 Kto. 1387 140831
 IBAN: DE56 2005 0550 1387 1408 31
 SWIFT / BIC HASPDE33XXX

