

Technical Datasheet

		EN ISO 13485 certified
Article - No:	11-222-0000	Revision: 2
Remarks		single use only

Specification

Product	Mailing container PP, 50ml, Ø44x114 mm, with absorbing inlay and screw thread, without cap, natural, autocl. up to 121°C
Dimensions	diameter: 44,00 mm (+/- 0,2 mm)
	height: 114,00 mm (+/- 0,2 mm)
Colour	natural
Nominal Capacity	50ml
Material Characteristics	container: polypropylene (PP)
	inlay: cellulose
Temperature range	-196°C - 121°C
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Packaging	pcs per bag: 250
	pcs per carton: 250
	pcs per pallet: 5.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

