

Technical Datasheet

		EN ISO 13485 certified
Article - No:	11-243-1200	Revision: 4
Remarks		single use only

Specification

Product	Mailing container PE, 30x30x184 mm, with enclosed screw cap PE, for swabs, natural
Colour	natural
Material Characteristics	container: high density polyethylene (HDPE) cap: polyethylene (PE)
Temperature range	-50°C – 70°C
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62
Packaging	pcs per bag: 375 pcs per carton: 375 pcs per pallet: 6.000
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C pal. not stackable
Approved by nerbe plus quality department	