

Technical Datasheet

		EN ISO 13485 certified
Article - No:	11-343-0105	Revision: 2
Remarks		single use only

Specification

Product	Mailing bag LDPE, 170x310 mm, triple punched in two rows, opaque, dark blue
Dimensions	width: 170,00 mm (+/- 0,2 mm)
	length: 310,00 mm (+/- 0,2 mm)
Colour	dark blue
Material Characteristics	low density polyethylene (LDPE)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Packaging	pcs per bundle: 100
	pcs per carton: 1.000
	pcs per pallet: 48.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

nerbe plus GmbH & Co. KG
 Porschestr. 25
 D-21423 Winsen/Luhe
 Tel: +49 (0)4171 7849 0
 Fax: +49 (0)4171 6 45 45

info@nerbe-plus.de
 www.nerbe-plus.de
 HRA Lüneburg 203240
 USt.-ID: DE 32583 0652

Pers. haftende Gesellschafterin:
 NERBE R&D GmbH, Winsen/Luhe
 GF: Philipp J.F. Nerbe
 HRB Lüneburg 208826

Deutsche Bank AG
 Kto. 0204 41200
 IBAN: DE41 2407 0075 0020
 4412 00
 SWIFT / BIC DEUTDE2H240

Hamburger Sparkasse
 Kto. 1387 140831
 IBAN: DE56 2005 0550 1387 1408 31
 SWIFT / BIC HASPDEHHXXX

