

Technical Datasheet

EN ISO 13485 certified		
Article - No:	11-369-0400	Revision: 1
Remarks		single use only

Specification

Product	Mailing bag PAP, for diagnostic specimen, 212x132x33 mm, UN approved 3373, suitable for nerbe mailing container, incl. assembly instruction, brown
Dimensions	length: 212,00 mm (+/- 0,2 mm)
	width: 132,00 mm (+/- 0,2 mm)
	height: 33,00 mm (+/- 0,2 mm)
Colour	brown
Material Characteristics	cardboard (PAP)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Packaging	pcs per carton: 150
	pcs per pallet: 4.400
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	