

Technical Datasheet

		EN ISO 13485 certified
Article - No:	11-632-0000	Revision: 2
Remarks		single use only

Specification

Product	Screw cap container PP, 35ml, Ø39x55 mm, with enclosed screw cap PE, container transparent, cap natural
Dimensions	diameter: 39,00 mm (+/- 0,2 mm)
	height: 55,00 mm (+/- 0,2 mm)
Colour	container transparent, cap natural
Nominal Capacity	35ml
Material Characteristics	container: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	-10°C - 40°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	flexible and unbreakable
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 500
	pcs per carton: 500
	pcs per pallet: 12.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

00-000-0001/06.2020

