

Technical Datasheet

Date: 22.06.20

Page 1 of 2

EN ISO 13485 certified	
Article - No:	11-702-6110
Remarks	Revision: 1 single use only

Specification

Product	Screw cap container PP, 500ml, Ø90x91 mm, with fitted screw cap PP, transparent, sterile R
Dimensions	container:
	diameter: 89,20 mm (+/- 0,2 mm)
	height: 90,50 mm (+/- 0,2 mm)
	cap:
	diameter: 96,00 mm (+/- 0,2 mm)
	height: 15,70 mm (+/- 0,2 mm)
	total height container and cap: 91,00 mm (+/- 0,2 mm)
Colour	transparent
Nominal Capacity	500ml
Material Characteristics	polypropylene (PP)
Temperature range	-80°C - 70°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
	flexible and soft
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 160
	pcs per carton: 160
	pcs per pallet: 1.280
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C

00-000-0001/06.2020



Technical Datasheet

Date: 22.06.20

Page 2 of 2

pal. not stackable

Approved by nerbe plus quality department
--