

Technical Datasheet

		EN ISO 13485 certified
Article - No:	11-822-6100	Revision: 2
Remarks		single use only

Specification

Product	Multi purpose container PP, 400ml, Ø88x100 mm, with fitted snap-on lid PE, transparent, graduated, sterile R
Dimensions	outer diameter: 88,00 mm (+/- 0,2 mm)
	inner diameter: 84,40 mm (+/- 0,2 mm)
	height (with lid): 101,00 mm (+/- 0,2 mm)
	height (without lid): 100,40 mm (+/- 0,2 mm)
Colour	transparent
Graduation	grad.
Nominal Capacity	400ml
Material Characteristics	container: polypropylene (PP)
	cap: polyethylene (PE)
Temperature range	-80°C - 75°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	do not contain any pyrogens
Shelf - Life	36 months at controlled storage conditions
Packaging	pcs per bag: 185
	pcs per carton: 185
	pcs per pallet: 2.220
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

