

Technical Datasheet

		EN ISO 13485 certified
Article - No:	11-881-0052	Revision: 4
Remarks		single use only

Specification

Product	Drosophila vial PS, 50ml, Ø28,5x95 mm, Flachboden 100/Tray, highly transparent
Dimensions	Länge: 30 cm
	Breite 31 cm
	Höhe 10 cm
Colour	highly transparent
Nominal Capacity	50ml
Material Characteristics	polystyrene (PS)
Temperature range	-10°C – 70°C
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 100
	pcs per carton: 500
	pcs per pallet: 15.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	