

Technical Datasheet

Date: 22.06.20

Page 1 of 1

		EN ISO 13485 certified
Article - No:	11-881-1001	Revision: 2
Remarks		single use only

Specification

Product	Foam stopper for drosophila vial Ø 25 and 29 mm, Ø30x25 mm, mite-proof
Dimensions	diameter: 30,00 mm (+/- 0,2 mm)
	height: 25,00 mm (+/- 0,2 mm)
Material Characteristics	polyurethane (PUR)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
Packaging	pcs per cs: 2.500
	pcs per carton: 2.500
	pcs per pallet: 50.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	