

Technical Datasheet

		EN ISO 13485 certified
Article - No:	11-952-6000	Revision: 1
Remarks	filling level line at 33 ml	single use only

Specification

Product	Flip top container PP, 50ml, Ø29x85 mm, flat bottom, natural, np pcr ready, sterile R, autoclavable up to 121°C
Dimensions	outer diameter (bottom): 29,00 mm (+/- 0,2 mm)
	outer diameter (top): 31,20 mm (+/- 0,2 mm)
	inner diameter (top): 29,00 mm (+/- 0,2 mm)
	outer diameter (lid): 35,00 mm (+/- 0,2 mm)
	height (without lid): 83,00 mm (+/- 0,2 mm)
	height (with lid): 85,00 mm (+/- 0,2 mm)
Colour	natural
Nominal Capacity	50ml
Material Characteristics	polypropylene (PP)
Temperature range	-40°C - 70°C
Autoclavability	autocl. up to 121°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	24 months at controlled storage conditions
Packaging	pcs per bag: 600
	pcs per carton: 600
	pcs per pallet: 14.400
Storage	dry condition: 40%-60% rel. humidity



Technical Datasheet

Date: 22.06.20

Page 2 of 2

	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	