

Technical Datasheet

EN ISO 13485 certified		
Article - No:	11-982-6000	Revision: 3
Remarks	filling level line at 161 ml	single use only

Specification

Product	Flip top container PP, 270ml, Ø65x76 mm, flat bottom, natural, graduated, np pcr ready, sterile R, autoclavable up to 121°C
Dimensions	outer diameter (bottom): 65,00 mm (+/- 0,2 mm)
	outer diameter (top): 78,00 (+/- 0,2 mm)
	inner diameter (top): 74,7 mm (+/- 0,2 mm)
	outer diameter (lid): 81,50 mm (+/- 0,2 mm)
	height (without lid): 74,70 mm (+/- 0,2 mm)
	height (with lid): 76,00 mm (+/- 0,2 mm)
Colour	natural
Graduation	grad.
Nominal Capacity	270ml
Material Characteristics	polypropylene (PP)
Temperature range	-40°C - 70°C
Autoclavability	autocl. up to 121°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	free of PCR inhibitors
	free of human DNA
	free of RNase and DNase
	free of pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	24 months at controlled storage conditions
Packaging	pcs per bag: 120
	pcs per carton: 120
	pcs per pallet: 2.400

00-000-0001/06.2020



Technical Datasheet

Date: 22.06.20

Page 2 of 2

Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	