

Technical Datasheet

EN ISO 13485 certified		
Article - No:	12-001-5002	Revision: 4
Remarks	Universal fit	single use only

Specification

Product	GTS-pipette, Single-channel, variable, 1.000-10.000µl, color code turquoise, CE
Classification	CE
Colour	color code turquoise
Nominal Capacity	1.000-10.000µl
Purity Grade	manufactured in clean rooms of ISO Class 8 or better (DIN EN ISO 14644-1)
Shelf - Life	not limited at controlled storage conditions
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	