

Technical Datasheet

		EN ISO 13485 certified
Article - No:	12-113-0055	Revision: 5
Remarks	Graduation at 0,5ml / 1ml / 1,5ml/ 2ml / 2,5ml / 3ml (for guidance only)	single use only

Specification

Product	Transfer pipette PE, 3ml, length 160 mm, transparent, grad.
Dimensions	diameter (stem): 8,00 mm (+/- 0,5 mm)
	outer diameter (orifice): 3,00 mm (+/- 0,2 mm)
	diameter (bulb): 13,00 mm (+/- 0,5 mm)
	length (stem): 118,00 mm (+/- 0,5 mm)
	length (bulb): 42,00 mm (+/- 0,5 mm)
	total length: 160,00 mm (+/- 10 mm)
Colour	transparent
Graduation	grad.
Nominal Capacity	3ml
Material Characteristics	low density polyethylene (LDPE)
Temperature range	-50°C - 55°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per inner cs: 500
	pcs per carton: 3.000
	pcs per pallet: 90.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	

