

Technical Datasheet

Date: 25.09.23

Page 1 of 2

		EN ISO 13485 certified
Article - No:	12-113-8100	Revision: 3
Remarks	Graduation at 0,5ml / 1ml / 1,5ml/ 2ml / 2,5ml / 3ml (for guidance only)	single use only

Specification

Product	Transfer pipette PE, 3ml, single packed, length 150 mm, transparent, grad., sterile R
Dimensions	diameter (stem): 8,00 mm (+/- 0,5 mm)
	outer diameter (orifice): 3,00 mm (+/- 0,2 mm)
	diameter (bulb): 13,00 mm (+/- 0,5 mm)
	length (stem): 109,00 mm (+/- 0,5 mm)
	length (bulb): 41,00 mm (+/- 0,5 mm)
	total length: 150,00 mm (+/- 10 mm)
Colour	transparent
Graduation	grad.
Nominal Capacity	3ml
Material Characteristics	low density polyethylene (LDPE)
Temperature range	-50°C - 55°C
Product Characteristics	biol. inert
	physiol. harmless
	no heavy metals acc. to EC/94/62
	Food contact: raw material in compl. with reg. EC 1935/2004
Purity Grade	do not contain any pyrogens
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per bag: 200
	pcs per carton: 1.000
	pcs per pallet: 30.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable



Technical Datasheet

Date: 25.09.23

Page 2 of 2

Approved by nerbe plus quality department