

# Technical Datasheet

		EN ISO 13485 certified
<b>Article - No:</b>	<b>12-531-6106</b>	<b>Revision: 1</b>
<b>Remarks</b>		<b>single use only</b>

## Specification

<b>Product</b>	Spreader PS, Ø3,5x145 mm, width 38 mm, L-shaped, 5 pcs/bag, white, sterile R
<b>Dimensions</b>	diameter: 3,50 mm (+/- 0,2mm)
	length: 145,00 mm (+/- 0,2mm)
	width spatula: 38,00 mm (+/- 0,2mm)
<b>Colour</b>	white
<b>Material Characteristics</b>	polystyrene (PS)
<b>Product Characteristics</b>	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
	Food contact: raw material in compliance with regulation EC 1935/2004
<b>Purity Grade</b>	do not contain any pyrogens
<b>Sterilisation Process</b>	radiation sterilization SAL 10-6
<b>Shelf - Life</b>	36 months at controlled storage conditions
<b>Packaging</b>	pcs per bag: 400
	pcs per carton: 4.000
	pcs per pallet: 128.000
<b>Storage</b>	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
<b>Approved by nerbe plus quality department</b>	

00-000-0001/06.2020

