

Labels

Name	_____
Date	_____ nerbe plus

02-918-3006

Name	_____
Date	_____ nerbe plus

02-938-3006

Name	_____
Date	_____ nerbe plus

02-958-3006

All tubes can be provided with the appropriate labels from the factory for an additional charge. These labels can be delivered with printed-on additional information (e. g. the institution's address) for an additional charge. Other labels with a customized layout are available on request.

24-hour Urine Container	
<ol style="list-style-type: none"> 1. Empty bladder into the toilet and note down time. 2. Then collect a sample of urine for each urination throughout the day and the following night into the quadrangular urine container. Empty urine into the container before each bowel movement. 3. On the next morning empty urine into the container at the same time as noted the day before. 4. Mix urine well. 5. Measure the amount of urine at the scale and write the amount on the labeling field of the urine container or on a smaller transport box. 6. If necessary, fill a part of the urine into the transporter box and give to the physician. 7. Store urine container in the refrigerator. Bottle and lid close tightly and do not let any urinary odor escape. 	
Store in a cool place!	
Date	_____
Name	_____
DOB	_____
Ward	_____
Amount	_____ ml
from	_____ to _____ (time)
nerbe plus	
D-21423 Winsen/Luhe, Germany	

02-978-3006

Order no.	Size in mm	Ground color	Impression	Packaging unit carton	Remarks
02-918-3006	50 x 25	white	black	10,000	Minimum purchase: 10,000 pcs.
02-938-3006	65 x 35	white	black	10,000	Minimum purchase: 10,000 pcs.
02-958-3006	75 x 20	white	black	10,000	Minimum purchase: 10,000 pcs.
02-978-3006	76 x 168	white	black		with instructions for the collection of 24-hour urine
02-998-3006	100 x 140	white	purple	1,000	Labels for envelopes with printed-on Aesculapian staff. With printed-on address on request.
02-902-3007	50.8 x 5.6	white	black		Barcode label for Thermoquick plate DMSO-resistant, freeze-resistant up to -80 °C