

**np** nerbe plus

GERMANY

---

SEROLOGISCHE PIPETTEN / SEROLOGICAL PIPETTES



ICONS  
ICONS



Artikelnummer  
Order no.



Information  
Information



Material  
Material



Neues Produkt  
New product



Innenkarton  
Box



Abmessungen  
Measurements



Umkarton  
Master shipper



Steril R (strahlensteril)  
Sterile R (sterile by radiation)



Polybeutel  
Bag



**np** herbe plus

10 ml in 1/10  
STERILE IF UNOPENED  
UNDAMAGED

10ml in 1/10  
TPLEX20CSEROTOCAL

herbe plus  
10ml in 1/10

herbe plus  
10ml in 1/10

herbe plus  
10ml in 1/10

herbe plus  
10ml in 1/10

herbe plus  
10ml in 1/10

# SEROLOGISCHE PIPETTEN

## SEROLOGICAL PIPETTES

- **Fussel freier Filter aus PE**  
Lint-free PE filter
- **Nicht zytotoxisch**  
Non-cytotoxic
- **Farbcodiert**  
Color coded




REF	MAT						
12-521-9101	PS	275 mm	2 ml	50		500	Aspirationspipette, Einzel-Peelpack, ohne Farbcode / aspirating pipette, single peelpack, no color code <b>STERILE R</b>
12-321-6103	PS	270 mm	1,1 ml	25		1.000	Farbcode: gelb / color code: yellow <b>STERILE R</b>
12-331-6103	PS	275 mm	2,2 ml	25		1.000	Farbcode: grün / color code: green <b>STERILE R</b>
12-401-6103	PS	270 mm	1 ml	25		1.000	Farbcode: gelb / color code: yellow <b>STERILE R</b>
12-401-9103	PS	270 mm	1 ml			500	Einzel-Peelpack, Farbcode: gelb / single peelpack, color code: yellow, Dispenser-Box <b>STERILE R</b>
12-421-6101	PS	275 mm	2 ml	25		1.000	Farbcode: grün / color code: green <b>STERILE R</b>
12-421-9101	PS	275 mm	2 ml			500	Einzel-Peelpack, Farbcode: grün / single peelpack, color code: green, Dispenser-Box <b>STERILE R</b>



# SEROLOGISCHE PIPETTEN SEROLOGICAL PIPETTES

- **Fussel freier Filter aus PE**  
Lint-free PE filter
- **Nicht zytotoxisch**  
Non-cytotoxic
- **Farb codiert**  
Color coded




REF	MAT						
12-441-6105	PS	335 mm	5 ml	25		500	Farbcode: blau / color code: blue <b>STERILE R</b>
12-441-9105	PS	335 mm	5 ml			200	Einzel-Peelpack, Farbcode: blau / single peelpack, color code: blue, Dispenser-Box <b>STERILE R</b>
+ 12-441-9115	PS	335 mm	5 ml	25		300	Einzel-Peelpack, Farbcode: blau / single peelpack, color code: blue <b>STERILE R</b>
12-461-6108	PS	340 mm	10 ml	25		500	Farbcode: orange / color code: orange <b>STERILE R</b>
12-461-9108	PS	340 mm	10 ml			200	Einzel-Peelpack, Farbcode: orange / single peelpack, color code: orange <b>STERILE R</b>
+ 12-461-9118	PS	340 mm	10 ml	25		300	Einzel-Peelpack, Farbcode: orange / single peelpack, color code: orange, Dispenser-Box <b>STERILE R</b>
12-481-9102	PS	345 mm	25 ml			150	Einzel-Peelpack, Farbcode: rot / single peelpack, color code: red, Dispenser-Box <b>STERILE R</b>
+ 12-481-9112	PS	345 mm	25 ml	25		200	Einzel-Peelpack, Farbcode: rot / single peelpack, color code: red <b>STERILE R</b>
12-491-9102	PS	345 mm	50 ml			100	Einzel-Peelpack, Farbcode: lila / single peelpack, color code: violet, Dispenser-Box <b>STERILE R</b>

# SEROLOGISCHE PIPETTEN SPEZIALMODELLE


## SEROLOGICAL PIPETTES SPECIAL




12-441-7105




12-461-7108



12-471-8108  
12-471-9108



12-481-7105



REF	MAT						
12-441-7105	PS	230 mm	5 ml	50		200	Einzel-Peelpack, <b>kurze Ausführung</b> , Farbcode: blau / single peelpack, <b>short version</b> , color code: blue <b>STERILE R</b>
12-461-7108	PS	235 mm	10 ml	25		200	Einzel-Peelpack, <b>kurze Ausführung</b> , Farbcode: orange / single peelpack, <b>short version</b> , color code: orange <b>STERILE R</b>
12-471-8108	PS	340 mm	10 ml	25		500	<b>große Öffnung</b> , Farbcode: orange / <b>big orifice</b> , color code: orange <b>STERILE R</b>
12-471-9108	PS	340 mm	10 ml			200	Einzel-Peelpack, <b>große Öffnung</b> , Farbcode: orange / single peelpack, <b>big orifice</b> , color code: orange <b>STERILE R</b>
12-481-7105	PS	230 mm	25 ml	25		100	Einzel-Peelpack, <b>kurze Ausführung</b> , Farbcode: rot / single peelpack, <b>short version</b> , color code: red <b>STERILE R</b>



**nerbe plus GmbH & Co.KG**  
**Porschestraße 25**  
**21423 Winsen · Germany**  
**Tel: +49 (0)4171 7849 0**  
**Fax: +49 (0)4171 645 45**  
**E-Mail: [info@nerbe-plus.de](mailto:info@nerbe-plus.de)**  
**Web: [www.nerbe-plus.de](http://www.nerbe-plus.de)**