

Technical Datasheet

Date: 15.05.18

Page 1 of 1

		EN ISO certified
Article - No:	02-623-0000	Revision: 3
Remarks		single use only

Specification

Product	Pressure plug PE, Ø11 mm, natural
Dimensions	diameter: 11,00 mm (+/- 0,2 mm)
Colour	natural
Material Characteristics	polyethylene (PE)
Product Characteristics	biologically inert
	physiologically harmless
	no heavy metals according to EC/94/62
Shelf - Life	not limited at controlled storage conditions
Packaging	pcs per bag: 1.000
	pcs per carton: 3.000
	pcs per pallet: 640.000
Storage	dry condition: 40%-60% rel. humidity
	bags / cases orig. closed
	storage temp. +10°C to max. +30°C
	pal. not stackable
Approved by nerbe plus quality department	