

Technical Datasheet

		EN ISO 13485 certified
Article - No:	06-377-2018	Revision: 10
Remarks		single use only

Specification

Product	Pipette tips PP, premium surface, refill system, 0-300µl, extra long, 96 pcs/tray, transparent tips, translucent yellow tray, grad., np premium bio, sterile R, autocl. up to 121°C, CE/IVD
Classification	CE/IVD
Dimensions	length: 59,35 mm (+/- 0,2 mm) inner diameter: 0,52 mm (+/- 0,2 mm)
Colour	transparent tips, translucent yellow tray
Graduation	grad.
Nominal Capacity	0-300µl
Material Characteristics	polypropylene (PP), premium surface
Autoclavability	autocl. up to 121°C
Product Characteristics	biologically inert physiologically harmless no heavy metals according to EC/94/62 sterile R SAL 10-6
Purity Grade	free of PCR inhibitors free of human DNA free of RNase and DNase free of ATP free of pyrogens and endotoxins
Sterilisation Process	radiation sterilization SAL 10-6
Shelf - Life	60 months at controlled storage conditions
Packaging	pcs per inner cs: 960 pcs per carton: 5.760 pcs per pallet: 345.600
Storage	dry condition: 40%-60% rel. humidity bags / cases orig. closed storage temp. +10°C to max. +30°C



Technical Datasheet

Date: 07.12.20

Page 2 of 2

pal. not stackable

Approved by nerbe plus quality department
--