


# np Safelock-Cap tube



The new optimal lid-mechanism of our **np Safelock-Cap** tube avoids accidental opening of the lid during centrifugation and thermal stress.

Sealing of the lid is guaranteed

- with **up to 30.000 x g** centrifugal impact,
- during thermal exposure in incubation or in a boiling water bath **up to 100° C**.

**All users were additionally convinced by the following benefits and advantages:**

### “HIGH SPEED PROTECTION”

Exceptional stability and highest safety in your centrifuge.

### “PRECISE LID-FIT”

Hardly any evaporation and sample loss thanks to high quality lid-fit.

### “EASY SAMPLE IDENTIFICATION”

Frosted writing area on the lid and tube.

### “CRYSTAL CLEAR”

Extremely clear product thanks to selected and particularly clear polypropylene, for perfect optical inspection of your sample material.

### “EASY TO USE”

Especially user-friendly thanks to our exceptionally designed lid rim.

### “AUTOCLAVABLE”

121° C, 15 psi, 15 minutes

Order no.	Volume	Material	Color	Packaging unit bag / carton
04-332-1200 / -1201 / -1202	0,5 ml	PP	clear / green / red	1.000 / 10.000
04-332-1203 / -1205 / -1206	0,5 ml	PP	yellow / blue / amber	1.000 / 10.000
04-212-1200 / -1201 / -1202	1,5 ml	PP	clear / green / red	1.000 / 8.000
04-212-1203 / -1205 / -1206	1,5 ml	PP	yellow / blue / amber	1.000 / 8.000
04-232-1200 / -1201 / -1202	2,0 ml	PP	clear / green / red	1.000 / 8.000
04-232-1203 / -1205 / -1206	2,0 ml	PP	yellow / blue / amber	1.000 / 8.000